ABOUT US

Xiamen Wiz Biotech Co., Ltd. is a high-tech biomedical company, devoted to the field of rapid diagnostic reagents and instrument. Wiz is located in Xiamen, China, established in July 2013, which is a listed enterprises in National Equities and Quotations. As a technology-driven company that prides itself on its scientific excellence, Wiz focused on technological innovation and product innovation with 15 authorized patents, each of our products embodies the creativity of our excellent scientists, who are working hard to continuously bring novel products to the China and world markets. SARS-CoV-2 rapid detection series of products, including SARS-CoV-2 antigen, antibody and influenza differential detection, are produced for Wizbiotech, which can easily and quickly carry out screening of SARS-CoV-2 in a large number of people, adding strength to epidemic prevention and control.

Self-testing

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)

- ✓ Anterior nasal swab specimen, noninvasive
- ✓ Very simple to use
- ✓ Convenient, no instrument required
- ✓ Rapid, results within 15~20 minutes
- ✓ Cost-efficient

<http://en.wizbiotech.com/>


Xiamen Wiz Biotech Co., Ltd.

Address: 3-4 Floor, NO.16 Building, Bio-medical Workshop, 2030 Wengjiao Xi Road,
Haicang District, Xiamen City, Fujian Province, 361026, P.R. China

Tel: +86(592)6808 278

E-mail: sales@wizbiotech.com

QUICK REFERENCE GUIDE

Video Operation

CLINICAL PERFORMANCE

Test Results	Reference PCR Results		
	Positive	Negative	Total
Positive	113	0	113
Negative	2	456	458
Total	115	456	571

Sensitivity: **98.26%** (95% C.I. 93.86%~99.79%)
 Specificity: **100.00%** (95% C.I. 99.19%~100.00%)
 Positive Predictive Value: **100%** (95% C.I. 96.79%~100.00%)
 Negativity Predictive Value: **99.56%** (95% C.I. 98.43%~99.95%)
 Overall Percent Agreement: **99.65%** (95% C.I. 98.74%~99.96%)

PACKING INFORMATION

Specification	Box size(mm)	Carton size(mm)	Carton G.W.(kg)	Boxes/Carton	Tests/Carton
1 Test	150*70*15mm	615*315*370mm	14kg	400	400
2 Tests	150*65*30mm	675*320*320mm	11kg	200	400
5 Tests	170*80*40mm	420*360*340mm	9kg	80	400
10 Tests	150*100*80mm	420*320*420mm	7kg	40	400
20 Tests	210*140*70mm	440*300*370mm	7kg	20	400

CERTIFICATE

CERTIFICATE

EC Certificate No. 1434-IVDD-492/2021

EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices

Polish Centre for Testing and Certification certifies
that manufactured by:

**Xiamen Wiz Biotech Co., Ltd.,
3-4 Floor, NO.16 Building, Bio-medical Workshop,
2030 Wengjiao Xi Road, Haicang District, Xiamen City,
Fujian Province, 361026, P.R. China**

in vitro diagnostic medical devices
for self-testing

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)


51332801, 51332802, 51332803, 51332804

in terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 22.11.2021 to 27.05.2024

The date of issue of the Certificate: 22.11.2021

The date of the first issue of the Certificate: 22.11.2021

Issued under the Contract No. MD-77/2021
Application No: 130/2021
Certificate bears the qualified signature.
Warsaw, 22/11/2021
Module A1

Vice-President

POLISH CENTRE FOR TESTING AND CERTIFICATION 02-844 Warsaw, 469 Pulawska Street, tel. +48 22 46 45 200, e-mail: pbc@pbc.gov.pl

FACTORY

Our one hundred thousand grade clean and dust-free workshop is operated strictly under ISO 13485:2016 and GMP guidelines. With the automatic production line, efficient production process, strict quality control, we always produce the products timely to meet customers' request.

INSTRUCTIONS FOR USE

EN

EN

SARS-CoV-2 Antigen Rapid Test(Colloidal Gold)

Instructions for Use

For use at physician's office or Non-professional
For use in non-clinical setting (e.g. home) nasal swab specimen

INTENDED USE

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) is intended for the qualitative detection of SARS-CoV-2 Antigen. Antigen-colloidal protein which is in nasal cavity (anterior nasal swab specimen from individuals with suspect COVID-19 infection). The test kit is intended for self-test or home test.

SUMMARY

Currently the main cause of COVID-19 infection is contact with someone who is already infected with the SARS-CoV-2, and asymptomatic infected people may also transmit the virus. Studies have shown that symptoms of infection generally appear within 14 days, with most occurring within 3 to 7 days after infection. The main symptoms are fever, fatigue, loss of smell and/or dry cough. In some cases, a stuffy nose, runny nose, muscle pain and diarrhea can also occur.

PRINCIPLE THE DETECTION

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) uses the colloidal gold method to detect the SARS-CoV-2 antigen in the test specimen. SARS-CoV-2, the SARS-CoV-2 Antigen will react with the antibody coated on the test line (T) to appear a red band in the test line (T) area; when the content of the SARS-CoV-2 in the test specimen is too low or does not exist, the test line (T) does not appear a red band. Regardless of whether the test specimen contains the SARS-CoV-2 antigen, there will be a red band in the area of the quality control line (C), which is the basis for judging whether the test is effective.

MAIN KIT COMPONENTS

Catalogue number	51332801	51332802	51332803	51332804
Specification	1 Testkit/bag	2 Testkit/bag	5 Testkit/bag	10 Testkit/bag
Test tube	1	2	5	10
Extraction tube	1	2	5	10
Disposable swab	1	2	5	10
Bioboard Waste Bag	1	2	5	10
Instructions for Use	1	1	1	1

Materials required but not provided: Timer

PRECAUTIONS

- For in Vitro Diagnostic Use Only.
- For use with nasal cavity (anterior nasal) swab specimen.
- For the detection of protein from SARS-CoV-2 only, not for detect other viruses or pathogens.
- When using this test kit, you must be aware of the law, the assistance of the legal guardian or authorized person.
- Keep the test kit away from children or pets.
- The test kit contains desiccant, it is forbidden to eat.
- The sample extraction solution in extraction tube contains chemical components. Direct contact should be avoided and eat prohibited. If the solution contacts the skin or mucous or eye, please flush with copious amounts of water. Please contact your family doctor or professional or seek medical advice if necessary.
- Use of gloves or other protective equipment are recommended when conducting testing.
- The test kits should be stored according to storage conditions required in Instruction for Use. It is recommended that the test kit is not stored in direct sunlight.
- Don't use the test kit beyond its expiration date.
- Don't use kit component that have been opened or changed.
- Never leave test card sealed in its foil pouch until just before use. Don't use if pouch is damaged or open.
- Do not use if the test kit or its foil pouch packaging is damaged or opened.
- The use of the swab should strictly follow instruction for Use, otherwise it may cause nasal cavity bleeding, swab rupture and retention and other risks.
- Don't immerse the swab in the extraction solution or other liquid before collect sample by disposable swabs.
- Don't touch swab soft tip when handling the swab sample.
- Proper sample collection and handling are crucial for correct results.
- Don't mix components from different kit lots.

TEST PROCEDURE

It is important to read the Instructions for Use carefully and follow the steps in the correct order. It takes about 15 minutes to prepare before each test, and the results can be obtained after 30 minutes.

Test Procedure	Preparatory work	Sample collecting	Sample processing	Sample testing
Chapter	I	II	III	IV

I Preparatory work

- Please use the test kit at room temperature (15°C~30°C). If the test kit was previously stored in a cold place (temperature below 15°C), please balance it at 15~30°C for 30 minutes before test.

- Prepare a Timer (such as watch, clock, tissue, hand sanitizer) and warm water.

- Please read the instruction for Use carefully.

- Wash hands thoroughly (at least 20 seconds) with soap and warm water. Hand sanitizer. This guarantees that the kit will not be contaminated, then dry your hands.

- Check the components in the box and verify that there is no damage or breakage.

Test device (sealed package)	Disposable swab (sealed package)	Extraction tube (filled with extraction solution)	Bioboard Waste Bag	Timer (prepare by yourself)

II Sample collecting

- Take out the sample extraction tube; unscrew the tube cap of the extraction tube.

- Tighten the tube cap for standby.

- Place the extraction tube on the bracket (attached to the box) to avoid liquid spillage.

- Tear off the swab's packing from stick end and take out the swab.

- As shown in the picture, wipe two nostrils with a swab.
 (1) Insert swab head into nostril less than 1 inch (Normally about 0.5~0.75 inch).
 (2) Gently rotate the swab against the nasal wall with moderate force at least five times.
 (3) Use the swab to repeat the collection procedure in another nostril.

- After two nostrils are collected, gently squeeze the extraction tube, and drop 2 drops liquid vertically into the sample well of the test card.

- Gently squeeze the extraction tube, and drop 2 drops liquid vertically into the sample well of the test card.

- Mix test kit is completed, put all test kit materials into the bioboard waste bag and dispose it according to the local bioboard waste disposal policy.

- Rinse hands thoroughly (at least 20 seconds) with soap and warm water/hand sanitizer.

IV Sample testing

- Tear off the aluminum foil bag, take out the test card and place it horizontally on the test desk.

⚠️ Warning!

The platform should be in a horizontal and stable state, and tilt and shake are strictly prohibited.

- Unplug the adding sample hole cover of the extraction tube.

⚠️ Warning!

The existence of bubbles in the extraction tube may lead to the wrong sample volume and inaccurate test result. If there are bubbles in the extraction tube, gently shake the extraction tube to squeeze out part of the liquid so as to remove the bubbles.

- Start timing, read the test results at 15 minutes. If the result is negative, read the test results at 30 minutes. Do not read result before 15 minutes or after 30 minutes.

- After test kit is completed, put all test kit materials into the bioboard waste bag and dispose it according to the local bioboard waste disposal policy.

- Rinse hands thoroughly (at least 20 seconds) with soap and warm water/hand sanitizer.

INTERPRETATION OF THE RESULTS

Negative:

The quality control line (C line) appears red band, while the test line (T line) doesn't appear red band.

Positive:

Both the quality control line (C line) and the test line (T line) appear red bands.

⚠️ Caution!

Regardless of the shades of color within the specified detection time, the result should be judged as positive result.

Invalids

Once the red line in the control line (C line) region doesn't appear which will be treated as invalid.

⚠️ Warning!

Invalid results indicate that you did not strictly follow the operation steps of instructions for Use. Please read the instruction for use carefully again and select a new kit for testing again.

PROCESSING OF TEST RESULTS

■ If the test result is positive:

- COVID-19 infection is present currently.

- Contact your doctor or local health department immediately.

- Contact your employer for specific requirement and protection guidelines in your area.

- Carry out PCR test for confirmation.

■ If the test result is negative:

- You need to continue to take measures related to contact with others and self-protection.

- Even if your test result is negative, it is also possible that there is an infection.

- If you still suspect you may be infected, you should repeat the test after 1-3 days. On account of the coronavirus can be accurately detected at every stages of infection.

■ If the test result is invalid:

- It may be caused by incorrect operation in the detection process.

- If the test result is still invalid, please contact your doctor or COVID-19 testing center.

QUALITY CONTROL

The quality control line is a key point of test kit and is used to control the procedure.

The quality control line appears, which indicates that the test has been performed correctly and the test kit has reacted.

SYMBOLS

Symbol	Used for	Symbol	Used for	Symbol	Used for
	General instructions		In vitro diagnostic medical device		Do not use
	In vitro Diagnostic Medical Device		Barcode		Information
	Shelf life -30°C~30°C		Barcode number		Lot code
	Authorized Representative in the European Union		Caution		Keep dry
	Do not use for medical purposes		Biological risks		In live devices
	Indicates a medical device		Indicates a medical device		Indicates a medical device
	Information label		Information label		Information label

Xiamen Wis Biotech Co., Ltd.
Address: 3/Floor, NO.16 Building, Bio-medical Workshop, 2030 Wengjiao XI Road, Haicang District, Xiamen City, Fujian Province, 361026, P.R. China
Tel: +86-592-6802878 Fax: +86-592-680279

Quray EC-REP BV
Pax 257, 2440 Geel, Belgium

CE 1434

SUPPLIERS OF DISPOSABLE STERILE SWAB

Zhejiang Gongdeng Medical Technology Co.,Ltd.

CE 0123 (According to Directive 93/42/EEC)

Address: No.10 Beiyuan Ave.,Huangyan,318002 Taizhou,Zhejiang,PEOPLE'S REPUBLIC OF CHINA

EC REP Name: Shanghai International Holding Corp.GmbH(Europe)

EC REP Address: Eiffestrasse 80,20337 Hamburg,Germany

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[2] Xu Chao, Li Ran, Analysis on the Risk Management of In Vitro Diagnostic Reagents[J]. China Medical Device Information, 2020, 26(13):10.

[3] Wu Yiqiang, Yang Rui, Rapid test of monoclonal gold with membrane as solid phase carrier[J]. Progress in Microbiology and Immunology, 2013, 31(1):74-76.

[4] Li Hongqi, Yang Rui, Rapid test of monoclonal gold with membrane as solid phase carrier[J]. Progress in Microbiology and Immunology, 2013, 31(1):74-76.

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