

# CLEANING AND DISINFECTION SOLUTION FOR MEDICAL DEVICES IN DENTAL ENVIRONMENTS

## DENTASEPT® SH PRO

### DISINFECTANT DETERGENT SOLUTION

Dentasept SH Pro is a ready-to-use cleaning and disinfection product for medical devices in a dental environment.

It combines broad spectrum efficacy, short contact times, excellent material compatibility and the properties of a gentle cleaner.

Dentasept SH Pro is suitable for use on a wide range of medical devices and can be used in practice as intended in your protocols, delivering the disinfection results and compliance level needed to restrict the spread of pathogens.

- Can be used on alcohol sensitive devices
- Good foaming properties
- Efficient formula without unnecessary ingredients (perfume and dye free)
- Cleaning and disinfection in one step, leaving the device visibly clean



# DENTASEPT® SH PRO

Disinfectant detergent solution, intended for cleaning and disinfection of non-invasive and nonsubmersible medical devices in dental environments

## PRODUCT CHARACTERISTICS


- Cleans and disinfects in one step
- Without alcohol, perfume and dye
- Broad spectrum antimicrobial and optimised contact time
- Compact foam dispenser to limit the formation of aerosols and volatile components (750 ml spray bottle)
- Formula ensures compatibility with a wide range of materials, including polymers (e.g. polycarbonate, plexiglas)
- Dense foam for better adhesion (750 ml spray bottle)
- pH 6

## INSTRUCTIONS FOR USE

Ready to use.



Spray the disinfectant detergent solution on a non-woven cloth and wipe.



Respect the indicated contact time for the required antimicrobial activity.



Rinsing is unnecessary except for areas intended to come into contact with skin or mucosa.



Please consult manufacturer's instruction for use of cleaning and disinfection of medical devices. For the frequency of use and cleaning of the application equipment, refer to the hygiene plan in place.

## MICROBIOLOGICAL PROPERTIES

CLAIM/CONDITIONS	CONTACT TIME
Bactericidal* - According to: EN 16615 (2 min), EN 13727 (1 min)	2 min
Yeasticidal* - According to: EN 13624, EN 16615	2 min
Tuberculocidal* - According to EN 14348	60 min
Active against enveloped viruses* - According to EN 14476	1 min
Active against Polyomavirus SV 40* - According to EN 14476	5 min
Active against Rotavirus* - According to EN 14476	30 sec
Limited spectrum virucidal* - According to EN 14476	30 min
Adenovirus	30 min
Norovirus (MNV)	15 min

\* dirty conditions

## ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Dentasept SH Pro	12 x 750 ml spray bottle	XXXXXX
Dentasept SH Pro	1 x 5 L canister	XXXXXX



D.M.D.  
1, Rue de l'Espoir  
59260 Lezennes  
France

Distributed by:  
**LABORATOIRES ANIOS**  
1, Rue de l'Espoir  
59260 Lezennes  
France

## APPLICATION AREA

Disinfectant detergent solution in dental environment, intended:

- Between two patients for cleaning and disinfection of non-invasive and nonsubmersible medical devices.
- Before the step of cleaning, disinfection and sterilization of medical devices (for example, rotary hand pieces, turbines, contra-angles), and in compliance with the indication of the manufacturer for wrapping in case of transport to the sterilization room and/or for wiping.

## INSTRUCTIONS FOR USE

Ready to use.

- 1) Spray the disinfectant detergent solution on a non-woven cloth and wipe.
- 2) Respect the indicated contact time for the required antimicrobial activity.
- 3) Rinsing is unnecessary except for the areas intended to come into contact with the skin or mucosa.
- 4) Please consult manufacturer's instruction for use for cleaning and disinfection of medical devices. For the frequency of use and cleaning of the application equipment, refer to the hygiene plan in place.

## INGREDIENTS

Didecylidimethylammonium chloride, excipients.

## CERTIFICATES

Medical Device. Full antimicrobial activity for Dentasept SH Pro is included in the scientific dossier, available on request.

## STORAGE

Storage: between +5°C and +25°C.

## MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.



www.anios.com