

Data sheet Royal Dent Putty Soft Fast			
Manuf. Date	12/09/2016	Batch number	248128
Test conditions	Test results	Ref. Parameters*	Warrant spec.**
Room Temperature (°C)	23,4	(23 ± 1) °C	/
Relative Humidity (%)	49,8%	(50 ± 10) %	/
Base Temperature (°C)	23	(23 ± 1) °C	/
Catalyst Temperature (°C)	23	(23 ± 1) °C	/
Requirements	Test results	Ref. Parameters*	Warrant spec.**
Total Working Time (sec)	73,6	70 – 100	65 – 100
Setting Time (sec)	151,15	135 – 170	120 – 160
Hardness 8' (Shore A)	63	57 – 63	58 – 65
<p>Badia Polesine, 13/09/2016</p> <p>QUALITY CONTROL DEPARTMENT</p> <p>This document is electronically released and signature is not necessary</p> <p>*Values refer to the date of product release. Test is made in conformity to our internal procedures described in our internal Quality Sheets: mixing ratio in weight, product handmixed and at humidity and temperature conditions mentioned in this Data Sheet.</p> <p>** Values refer to the features of the product until expiry date.</p>			