

SILICONES - CERTIFICATE OF ANALYSIS

B220002 ROYAL DENT LIGHT FAST 2x50ML

MANUFACTURED DATE: 2025-04-24		BATCH NUMBER: 0000469015		EXPIRY DATE: 2028-04-22	
Test parameter		Test Result	Specification	Test method	
Room temperature (°C)		23,4	[22,0-24,0]	2100-WI QC-001	
Relative Humidity (%RH)		47,0	[40,0-60,0]	2100-WI QC-001	
Base temperature (°C)		23,1	[22,0-24,0]	2100-AM QC-023	
Catalyst temperature (°C)		23,3	[22,0-24,0]	2100-AM QC-023	
Mixing time (s)		Compliant	Automatic		